

Pioneer in Imaging
www.med-works.in



**MED
WORKS**
www.med-works.in

High Frequency Mobile Digital Radiography System

AIRBILE 100



High Frequency Mobile
Digital Radiography System

Korean Technology



AIRBILE 100

High Frequency Mobile /
Digital Radiography System



AIRBILE-100 The Advent of an Advanced Mobile System and Easy Handling Compact and Lightweight design to achieve easy manoeuvrability and good visibility



Simple and Precise Positioning

- The effortless arm positioning with free stop function
- Robust structure

Effective PC Cradle for DR Users

Excellent Image Quality

- Only 1.5mm tube focal spot, superior in sharpness than any others

Superb Battery Technology

- The ultralight system under 98kg by applying a patented battery
- High-density small battery type
- Full charge within only 3.5hours

Convenient Movement

- The sleek and compact design only 52cm width
- The great pivotability with 360° rotating wheel
- A convenient and safe footbrake
- DR supportable using bluetooth

Smart Operation

- Applying remote control
- Dual operation panels (main & tube head console) for more efficient workflow



Technical Specifications



5kW, 100mA High Frequency Generator



The Ultra-Light System
• Under 98kg
• Outstanding battery technology



Free Stop Arm



Slim & Compact Design
• Only 52cm width



1.5mm Tube Focal Spot



Max. Anode Heat Capacity
• 3kJ



Convenient Movement



Simple and Precise Positioning



Smart Operation



Superb Battery Technology



MEDWORKS SOLUTIONS PVT LTD

89, Mahia Industrial Estate, Uranganpatti,
Madurai - 625 020. Tamilnadu, India.

☎ 0452 - 2323307, 0452 - 2424077

☎ 93445 33077, 98421 33077

✉ Info@med-works.in

🌐 www.med-works.in

