

Pioneer in Imaging
www.med-works.in



**MED
WORKS**
www.med-works.in

High Frequency Mobile Digital Radiography System

AIRPOLE 300



High Frequency Mobile
Digital Radiography System

Korean Technology



Technical Specifications

Generator	Input Power	Single Phase 200 - 240 VAC
	Capacity	15KW
	kV Range	40-125kV
	mA Range	10-300mA
	mAs Range	0.1-250mAs
	Time Range	0.001~ 10 sec
	Frequency	40kHz
	kV Accuracy	± 10 %
	mA Accuracy	± 15 %
	Anatomical Progamming	216APR
Tube	Anode heat storage capacity	140KHU
	Focal Spots	1.0/2.0mm Dual
	Target Angle	16 °
	X-ray Tube Anode	Rotating
	Max Tube Voltage	125kV
Collimator	Max Tube Voltage	150kV
	Min.Inherent Filtration	1.2mmAL
Display	Digital Display	7 segment LED
Storage Space	Laptop, DR panel	



MEDWORKS SOLUTIONS PVT LTD

89, Mahia Industrial Estate, Uranganpatti,
Madurai - 625 020. Tamilnadu, India.

☎ 0452 - 2323307, 0452 - 2424077

☎ 93445 33077, 98421 33077

✉ Info@med-works.in

🌐 www.med-works.in

