

Pioneer in Imaging
www.med-works.in



**MED
WORKS**
www.med-works.in

Ultralight Portable / Digital X-Ray System

AIRRAY 20



Ultralight Portable / Digital X-Ray System

Korean Technology



Technical Specifications



1.6kW, 20mA, 1.0mm Focal spot

High Frequency Generator



Remote Controller for Exposure



Ultra-Light Portable X-ray Generator

5.3kgs



Hybrid Battery-powered Portable X-ray Generator

Wireless Charger for Battery in carrying case



Bluetooth USB Adapter

For F/W Download, PC (DR) Interlock



Dual operation panel

Main OP on Rear, Mini OP on Top



Ultralight Weight



Simple and Precise Positioning



Smart Operation



Superb Battery Technology



MEDWORKS SOLUTIONS PVT LTD

89, Mahia Industrial Estate, Uranganpatti,
Madurai - 625 020. Tamilnadu, India.

☎ 0452 - 2323307, 0452 - 2424077

☎ 93445 33077, 98421 33077

✉ Info@med-works.in

🌐 www.med-works.in

