

Pioneer in Imaging



**MED
WORKS**

www.med-works.in

High Frequency Portable Digital Radiography System

AIRWAY 40



High Frequency Portable
Digital Radiography System

Korean Technology



kV Range	40 ~ 100kV
Maximum mA	35mA
Maximum mAs	0.4 ~ 100mAs in 35 steps
Peak Power	2.4kW
Focal Spot Size	1.2 - 1.2 mm
Target Angle	16°
Technology	High Frequency
Anode Heat Storage	20kHu
Weight	14kg
Battery	14.8V, 5000mAh
Input Power	180-240V AC ,5A

- **100kV, 35mA, 2.4kW
1.2mm Focal Spot**
High Frequency Generator
- **Compact Portable X-ray
Generator**
14.0 Kgs
- **Display reverse Switch &
SID Measuring tape**
- **Excellence Battery Backup**



MEDWORKS SOLUTIONS PVT LTD

89, Mahia Industrial Estate, Uranganpatti,
Madurai - 625 020. Tamilnadu, India.

☎ 0452 - 2323307, 0452 - 2424077

☎ 93445 33077, 98421 33077

✉ Info@med-works.in

🌐 www.med-works.in

