

Pioneer in Imaging
www.med-works.in



**MED
WORKS**

www.med-works.in

Ultra Light Portable Dental Digital X-Ray Systems

DENT X



Ultra Light Portable Dental Digital X-Ray Systems

Korean Technology



kV Range	70kV
Maximum mA	2mA
Maximum mAs	2mAs
Focal Spot Size	0.4 x 0.4 mm
Target Angle	12.5°
Peak Power Output	140W
Anode Heat Storage	4300J
Display	OLED Display with soft switch
Weight	1.6kg
Battery	1450mAh/300 exposure per charge
Input Power	19VDC , 3.16A

- 70kV, 2mA, 140W,
0.4mm Focal Spot
- Compact & Light one hand
Operation to make exposure
- Ergonomical design reduces
risk to exposure by scattered
radiation
- Only 1.6kg with Battery
Providing upto 300 Exposure



MEDWORKS SOLUTIONS PVT LTD

89, Mahia Industrial Estate, Uranganpatti,
Madurai - 625 020. Tamilnadu, India.

☎ 0452 - 2323307, 0452 - 2424077

☎ 93445 33077, 98421 33077

✉ Info@med-works.in

🌐 www.med-works.in

