

Pioneer in Imaging
www.med-works.in



**MED
WORKS**
www.med-works.in

Ultra-Light Portable Dental X-Ray System

REXTAR X




Ultra-Light Portable Dental X-Ray System

Korean Technology



Technical Specifications

	Model	REXTAR X	
Input power	X-ray Unit	19 VDC 3.16A	
	Battery Charger	100-240 VAC, 50-60Hz	1.7A
Output power		140W	
Power Supply		11.1 VDC	
Battery		11.1 V , 1450 mAh	
Frequency		70kHz	
Tube kV/Ma		70kV/2mA(Fixed)	
Max. Deviation	icX	±4%	
	MAs	±10%	
X-ray Tube	Model Name	D-041	
	Inherent Filtration	Min. 1.0mmAl	
	Focal Spot size	0.4 X 0.4 mm	
Total filtration		1.5mm	
Weight		1.6kg	



MEDWORKS SOLUTIONS PVT LTD

89, Mahia Industrial Estate, Uranganpatti,
Madurai - 625 020. Tamilnadu, India.

☎ 0452 - 2323307, 0452 - 2424077

☎ 93445 33077, 98421 33077

✉ Info@med-works.in

🌐 www.med-works.in

