

Pioneer in Imaging
www.med-works.in



**MED
WORKS**


www.med-works.in

Ultra Light Portable Dental Digital X-Ray Systems

DENT Q



Ultra Light Portable Dental Digital X-Ray Systems



kV Range	60kV
Maximum mA	3mA
Maximum mAs	3mAs
Focal Spot Size	0.4mm IEC
Target Angle	19°
Peak Power Output	210W
Anode Heat Storage	4500J
Display	OLED Display with soft switch
Weight	< 2kg
Battery	2900mAh/ 250 exposure per charge
Input Power	22.2V DC

- 70kV, 3mA, 210W,
0.4mm Focal Spot
- Advanced High Frequency
type X-Ray System employing
Double Insulation Technology
- Compact Design providing
Easy handling
- Excellent Battery Technology
250 Exposure for one time
Charge



MEDWORKS SOLUTIONS PVT LTD

89, Mahia Industrial Estate, Uranganpatti,
Madurai - 625 020. Tamilnadu, India.

☎ 0452 - 2323307, 0452 - 2424077

☎ 93445 33077, 98421 33077

✉ Info@med-works.in

🌐 www.med-works.in

