

Pioneer in Imaging
www.med-works.in



**MED
WORKS**
www.med-works.in

Ultralight Portable / Digital X-Ray System

AIRWAY SERIES



Korean Technology



AIRWAY Mini



70 kV, 2mA, 140W,
0.4mm Focal Spot

High Frequency Generator



Ultra-Light Portable X-ray
Generator

1.6 Kg



High Quality Portable
X-ray System



Excellence Battery Backup

One time charge, provides
upto 300 Exposures



AIRWAY 10



90 kV, 10mA, 500W,
0.8mm Focal Spot

High Frequency Generator



Ultra-Light Portable X-ray
Generator

3.0 Kgs



Reduced Radiation Exposure
(AI Tech)

With over 50% reduction noise and scatter



Use Convenience

Patient Motion Alarm, Digital distance
sensor, Book mark

AIRRAY 20



1.6kW, 20mA, 1.0mm Focal spot

High Frequency Generator



Ultra-Light Portable X-ray Generator

5.0 Kgs



Bluetooth USB Adapter

For F/W Download, PC (DR) Interlock



Hybrid Battery-powered Portable X-ray Generator

Wireless Charger for Battery in carrying case



AIRWAY 40



**100kV, 35mA, 2.4kW
1.2mm Focal Spot**

High Frequency Generator



Compact Portable X-ray Generator

14.0 Kgs



**Display reverse Switch &
SID Measuring tape**



Excellence Battery Backup

Specifications	AIRWAY-MINI	AIRWAY-10
kV Range	70kV	50 ~ 90kV
Maximum mA	2mA	10mA
Maximum mAs	0.01~1.30 sec	0.1mAs - 10mA in 19 steps
Peak Power	140W	500W
Focal Spot Size	0.4x0.4 mm	0.4x0.8 mm
Target Angle	12.5°	16°
Technology	High Frequency	High Frequency
Anode Heat Storage	4300J	8kJ
Weight	1.6kg	3.0kg
Battery	11.1V, 1450mAh	22.2V, 2200mAh
Input Power	180-240V AC ,5A	180-240V AC ,5A

Specifications	AIRRAY-20	AIRWAY-40
kV Range	50 ~ 90kV	40 ~ 100kV
Maximum mA	20mA	35mA
Maximum mAs	0.4 ~ 20 mAs in 25 steps	0.4 ~ 100mAs in 35 steps
Peak Power	1.6kW	2.4kW
Focal Spot Size	1.0 mm	1.2 x 1.2 mm
Target Angle	16°	16°
Technology	High Frequency	High Frequency
Anode Heat Storage	7kJ	20kHu
Weight	5.0kg	14kg
Battery	3.7V, 5300mAh	14.8V, 5000mAh
Input Power	180-240V AC ,5A	180-240V AC ,5A



AIRWAY-MINI



AIRWAY-10



AIRRAY-20



AIRWAY-40



MEDWORKS SOLUTIONS PVT LTD

89, Mahia Industrial Estate, Uranganpatti,
Madurai - 625 020. Tamilnadu, India.

☎ 0452 - 2323307, 0452 - 2424077

☎ 93445 33077, 98421 33077

✉ Info@med-works.in

🌐 www.med-works.in

