

Pioneer in Imaging

www.med-works.in



**MED
WORKS**

www.med-works.in

High Frequency Portable Digital Radiography System

AIRWAY 10



High Frequency Portable Digital Radiography System

Korean Technology



Technical Specifications

kV Range	50 ~ 90kV
Maximum mA	10mA
Maximum mAs	0.1mAs - 10mAs in 19 steps
Focal Spot Size	0.8 mm
Target Angle	16°
Peak Power Output	500W
Anode Heat Storage	8kJ
Generator	High Frequency
Weight	3.0kg
Battery	22.2V, 2200mAh
Input Power	180-240V AC ,5A

- Reduced Radiation Exposure (AI Tech)
With over 50% reduction noise and scatter
- Use Convenience
Patient Motion Alarm, Digital distance sensor, Book mark
- Fast Diagnostics Solution (AI Tech)
9 major pulmonary diseases with over 94% accuracy within 3 sec
- Digital Distance indication
using Laser sensor



MEDWORKS SOLUTIONS PVT LTD

89, Mahia Industrial Estate, Uranganpatti,
Madurai - 625 020. Tamilnadu, India.

☎ 0452 - 2323307, 0452 - 2424077

☎ 93445 33077, 98421 33077

✉ Info@med-works.in

🌐 www.med-works.in

