

Pioneer in Imaging

[www.med-works.in](http://www.med-works.in)



**MED  
WORKS**

[www.med-works.in](http://www.med-works.in)

Ultra Light Portable Dental Digital X-Ray Systems

**DENT X**



Ultra Light Portable Dental Digital X-Ray Systems

Korean Technology



|                           |                                 |
|---------------------------|---------------------------------|
| <b>kV Range</b>           | 70kV                            |
| <b>Maximum mA</b>         | 2mA                             |
| <b>Maximum mAs</b>        | 2mAs                            |
| <b>Focal Spot Size</b>    | 0.4 x 0.4 mm                    |
| <b>Target Angle</b>       | 12.5°                           |
| <b>Peak Power Output</b>  | 140W                            |
| <b>Anode Heat Storage</b> | 4300J                           |
| <b>Display</b>            | OLED Display with soft switch   |
| <b>Weight</b>             | 1.6kg                           |
| <b>Battery</b>            | 1450mAh/300 exposure per charge |
| <b>Input Power</b>        | 19VDC , 3.16A                   |

- 70kV, 2mA, 140W,  
0.4mm Focal Spot
- Compact & Light one hand  
Operation to make exposure
- Ergonomical design reduces  
risk to exposure by scattered  
radiation
- Only 1.6kg with Battery  
Providing upto 300 Exposure



#### **MEDWORKS SOLUTIONS PVT LTD**

89, Mahia Industrial Estate, Uranganpatti,  
Madurai - 625 020. Tamilnadu, India.

☎ 0452 - 2323307, 0452 - 2424077

☎ 93445 33077, 98421 33077

✉ Info@med-works.in

🌐 www.med-works.in

