

Pioneer in Imaging



High Frequency Mobile / Digital Radiography System

# AIRBILE-100



## THE MOVEMENT AS LIGHT AS AIR

Korean Technology



# AIRBILE-100

High Frequency Mobile /  
Digital Radiography System



AIRBILE-100 The Advent of an Advanced Mobile System and Easy Handling Compact and Lightweight design to achieve easy manoeuvrability and good visibility

### Simple and Precise Positioning

- The effortless arm positioning with free stop function
- Robust structure

### Effective PC Cradle for DR Users



### Excellent Image Quality

- Only 1.5mm tube focal spot, superior in sharpness than any others

### Superb Battery Technology

- The ultralight system under 98kg by applying a patented battery
- High-density small battery type
- Full charge within only 3.5hours

### Smart Operation

- Applying remote control
- Dual operation panels (main & tube head console) for more efficient workflow

### Convenient Movement

- The sleek and compact design only 52cm width
- The great pivotability with 360° rotating wheel
- A convenient and safe footbrake
- DR supportable using bluetooth





**5kW, 100mA High Frequency Generator**



**Slim & Compact Design**  
• Only 52cm width



**The Ultra-Light System**  
• Under 98kg  
• Outstanding battery technology



**1.5mm Tube Focal Spot**



**Free Stop Arm**



**Max. Anode Heat Capacity**  
• 3kJ



Convenient Movement



Simple and Precise Positioning



Smart Operation



Superb Battery Technology



**MEDWORKS SOLUTIONS PVT LTD**

89, Mahia Industrial Estate, Uranganpatti,  
Madurai - 625 020. Tamilnadu, India.

☎ 0452 - 2323307, 0452 - 2424077

☎ 93445 33077, 98421 33077

✉ Info@med-works.in

🌐 www.med-works.in

