

Pioneer in Imaging



Ultralight Portable / Digital X-Ray System

AIRRAY-20



X - R A Y

THE MOVEMENT AS LIGHT AS AIR

Korean Technology




AIRRAY-20

Ultralight Portable /
Digital X-Ray System

- **Remote Controller for Exposure** (Option)


** 10 APR, DAP

** 1.5hr (Wire charging),
4.5hr (Wireless charging)



- **Wireless Charging PAD for Battery**

In carrying case



Charging

C Type magnetic connector

AIRRAY-20 The Advent of an Advanced Portable System and Easy Handling Compact and Lightweight design to achieve easy manoeuvrability and good visibility

Skin Guard

Easy attach / detachable



Bluetooth USB Adapter

F/W Download.
PC (DR) Interlock

Dual Operation Panel

* Main OP on Rear, Mini OP on Top





1.6kW, 20mA, 1.0mm Focal spot

High Frequency Generator



Remote Controller for Exposure



Ultra-Light Portable X-ray Generator

5.3kgs



Hybrid Battery-powered Portable X-ray Generator

Wireless Charger for Battery in carrying case



Bluetooth USB Adapter

For F/W Download, PC (DR) Interlock



Dual operation panel

Main OP on Rear, Mini OP on Top



Ultralight Weight



Simple and Precise Positioning



Smart Operation



Superb Battery Technology



MEDWORKS SOLUTIONS PVT LTD

89, Mahia Industrial Estate, Uranganpatti, Madurai - 625 020. Tamilnadu, India.

☎ 0452 - 2323307, 0452 - 2424077

☎ 93445 33077, 98421 33077

✉ Info@med-works.in

🌐 www.med-works.in

