


Ultralight Portable / Digital X-Ray System


# AIRRAY-20





**THE MOVEMENT AS LIGHT AS AIR**


 **1.6kW, 20mA, 1.0mm Focal spot**  
High Frequency Generator

 **Remote Controller for Exposure**

 **Ultra-Light Portable X-ray Generator**  
5.3kgs

 **Hybrid Battery-powered Portable X-ray Generator**  
Wireless Charger for Battery in carrying case

 **Bluetooth USB Adapter**  
For F/W Download, PC (DR) Interlock

 **Dual operation panel**  
Main OP on Rear, Mini OP on Top




# AIRRAY-20

Ultralight Portable /  
Digital X-Ray System

- **Remote Controller for Exposure**  
(Option)

- \*\* 10 APR, DAP
- \*\* 1.5hr (Wire charging),  
4.5hr (Wireless charging)




- **Wireless Charging PAD for Battery**

In carrying case



AIRRAY-20 The Advent of an Advanced Portable System and Easy Handling Compact and Lightweight design to achieve easy manoeuvrability and good visibility




### Charging

C Type magnetic connector

### Skin Guard

Easy attach / detachable



### Bluetooth USB Adapter

F/W Download.  
PC (DR) Interlock



### Dual Operation Panel

\* Main OP on Rear, Mini OP on Top

