

Pioneer in Imaging



Compact Portable Digital X-Ray System

PXP - 16 HF

PXP - 20 HF PLUS

PXP - 40 HF



Model	PXP-16HF	PXP-20HF Plus	PXP-40HF
Max. Output(kW)	1.2	1.6	2.4
kV Range per 1kV step	50-80	40-100	40-100
mA Range/No. of step	10, 15, 16	15, 20	16-35/5
mAs Range/No. of step	0.4-30/27	0.3-50/22	0.4-100/25
Tube Focal Spot(mm)	0.8 × 0.8	1.2 × 1.2	1.2 × 1.2
Target Angle	16°	16°	16°
Dual Laser Pointer	Standard	Optional	Optional
Power Supply	100-120V/200-240V	200-240V	100-120V/200-240V
No. of APR	6	4	8
Weight(kg)	5.3	9.8	12.8
Packing Dimension	51 x 36 x 37 (15 kg)	51 x 36 x 37 (19kg)	58 x 40 x 40 (21kg)



PXP - 16 HF



PXP - 20 HF PLUS



PXP - 40 HF



MEDWORKS SOLUTIONS PVT LTD

89, Mahia Industrial Estate, Uranganpatti,
Madurai - 625 020. Tamilnadu, India.

☎ 0452 - 2323307, 0452 - 2424077

☎ 93445 33077, 98421 33077

✉ Info@med-works.in

🌐 www.med-works.in

